



1. Identification

Product name	Clopidogrel Bisulfate
Product code	RLP004
CAS number	120202-66-6
Chemical name	Thieno [3,2-c] pyridine-5(4H)-acetic acid, alpha-(2-chlorophenyl)-6,7-dihydro-, methyl ester, (S)-, sulfate (1:1).
Suggested use	Specified quality tests and assay use only.
Advised limitations	Not for use as a drug. Not for administration to humans or animals.
Manufacturer	Reflabpharm LLC
Address	Armenia, Yerevan, 0019, Silikyan 2 nd str. 25/1
Telephone	+37411554055
Website	www.reflab.am
E-mail	info@reflab.am

2. Hazard(s) Identification

Physical hazards	Not classified.
Health hazards	Acute toxicity, oral Category 4 Skin corrosion/irritation Category 1 Serious eye damage/eye irritation Category 1 Specific target organ toxicity, single exposure Category 1 (Blood)
Environmental hazards	Not classified.
OSHA defined hazards	Not classified.
Label elements	
Signal word	Danger.
Hazard statement	Harmful if swallowed. Causes severe skin burns and eye damage. Causes damage to organs (Blood).
Precautionary statement	
Prevention	Do not breathe dust. Wash thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection. If exposed: Immediately call a poison center/doctor. If swallowed: Rinse mouth. Do NOT induce vomiting. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. Wash contaminated clothing before reuse. If inhaled: Remove person to fresh air and keep comfortable for breathing. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
Response	

Material name	Cas. N	Version #	Revision date	Issue date	Page
Clopidogrel Bisulfate	120202-66-6	2	20.05.2026	06.03.2024	1/8

Storage	Store locked up.
Disposal	Dispose of contents/container in accordance with local/regional/national/international regulations.
Hazard(s) not otherwise classified (HNOC)	None known.
Supplemental information	Pharmacologically active material.

3. Composition/information on ingredients

Substance				
Chemical name	Common name and synonyms	CAS number	%	
Clopidogrel Bisulfate	Thieno [3,2-c] pyridine-5(4H)-acetic acid, alpha-(2-chlorophenyl)-6,7-dihydro-, methyl ester, (S)-, sulfate (1:1).	120202-66-6	100	

4. First-aid Measures

Inhalation	If breathing is difficult, remove to fresh air and keep at rest in a position comfortable for breathing. Call a physician if symptoms develop or persist.
Skin contact	Take off immediately all contaminated clothing. Rinse skin with water/shower. Call a physician or poison control center immediately. For minor skin contact, avoid spreading material on unaffected skin.
Eye contact	Call a physician or poison control center immediately. Rinse cautiously with water for several minutes.
Ingestion	Rinse mouth. If ingestion of a large amount does occur, call a poison control center immediately. Do not induce vomiting without advice from poison control center. If vomiting occurs, keep head low so that stomach content doesn't get into the lungs.
Most important symptoms/effects, acute and delayed	Pharmacologically active material. Occupational exposure may cause physiological effects.
Indication of immediate medical attention and special treatment needed	Provide general supportive measures and treat symptomatically. Chemical burns: Flush with water immediately. While flushing, remove clothes which do not adhere to affected area. Call an ambulance. Continue flushing during transport to hospital.
General information	Remove from exposure. Remove contaminated clothing. For treatment advice, seek guidance from an occupational health physician or other licensed health-care provider familiar with workplace chemical exposures. If person is not breathing, give artificial respiration. If breathing is difficult, give oxygen if available. Persons developing serious hypersensitivity (anaphylactic) reactions must receive immediate medical attention.

5. Fire-fighting Measures

Material name	Cas. N	Version #	Revision date	Issue date	Page
Clopidogrel Bisulfate	120202-66-6	2	20.05.2026	06.03.2024	2/8

Suitable extinguishing media	Water. Foam. Dry chemical or CO ₂ . Use fire-extinguishing media appropriate for surrounding materials.
Unsuitable extinguishing media	None known.
Specific hazards arising from the chemical	No unusual fire or explosion hazards noted.
Special protective equipment and precautions for firefighters	Wear suitable protective equipment.
Fire-fighting instructions	Use water spray to cool unopened containers. As with all fires, evacuate personnel to a safe area.
Equipment	Firefighters should use self-contained breathing equipment and protective clothing.
Specific methods	Use standard firefighting procedures and consider the hazards of other involved materials.
General fire hazards	No unusual fire or explosion hazards noted.

6. Accidental Release Measures

Personal precautions, protective equipment and emergency procedures	Keep unnecessary personnel away. Do not touch damaged containers or spilled material unless wearing appropriate protective clothing. Ensure adequate ventilation. Wear appropriate personal protective equipment. Avoid inhalation of dust from the spilled material. For personal protection, see section "Personal Protection" of the SDS. Dust deposits should not be allowed to accumulate on surfaces, as these may form an explosive mixture if they are released into the atmosphere in sufficient concentration.
Methods and materials for containment and cleaning up	Avoid the generation of dusts during clean-up. Sweep up or vacuum up spillage and collect in suitable container for disposal. Clean surface thoroughly to remove residual contamination. For waste disposal, see section "Disposal considerations" of the SDS.
Environmental precautions	Avoid discharge into drains, water courses or onto the ground.

7. Handling and Storage

Precautions for safe handling	As a general rule, when handling CRMs, avoid all contact and inhalation of dust, mists, and/or vapors associated with the material. Clean equipment and work surfaces with suitable detergent or solvent after use. After removing gloves, wash hands and other exposed skin thoroughly. Select and use containment devices and personal protective equipment based on a risk assessment of material potency and exposure potential.
Conditions for safe storage, including any incompatibilities	This material should be handled and stored per label instructions to ensure product integrity.

8. Personal Protection

Occupational exposure limits	The following constituents are the only constituents of the product which have a PEL, TLV or other recommended
------------------------------	--

Material name	Cas. N	Version #	Revision date	Issue date	Page
Clopidogrel Bisulfate	120202-66-6	2	20.05.2026	06.03.2024	3/8

exposure limit. At this time, the other constituents have no known exposure limits.

Exposure limit values

Industrial Use Material	Type	Value
Clopidogrel Bisulfate (CAS 120202-66-6)	TWA	20 micrograms/m ³
Appropriate engineering controls	For laboratory operations, use local exhaust ventilation or a ventilated enclosure for high energy operations such as particle sizing. Control exposures to below the occupational exposure level (if available). Select and use containment devices and personal protective equipment based on a risk assessment of exposure potential. Cover all containers for solutions and slurries while being transferred.	
Biological limit values	No biological exposure limits noted for the ingredient(s).	

Individual protection measures

Eye/face protection	Wear safety glasses with side shields, chemical splash goggles, or full face shield, if necessary. Base the choice of protection on the job activity and potential for contact with eyes or face. An emergency eye wash station should be available.
Skin protection	
Hand protection	Wear nitrile or other impervious gloves if skin contact is possible. When the material is dissolved or suspended in an organic solvent, wear gloves that provide protection against the solvent.
Other	Train employees in proper gowning and degowning practices. Wear lab coat. Base the choice of skin protection on the job activity, potential for skin contact and solvents and reagents in use. Do not wear protective garments in common areas (e.g., cafeterias) or out-of-doors.
Respiratory protection	Respirators are generally not required for laboratory operations. Use a tight-fitting full-face respirator with HEPA filters for spill cleanup. Choose respiratory protection appropriate to the task and the level of existing engineering controls.
Thermal hazards	Wear appropriate thermal protective clothing, when necessary.
General hygiene considerations	Handling practices in this SDS are recommendations for laboratory use of reference standards. Procedures for any other uses or quantities should be determined after an appropriate assessment.

9. Physical and chemical properties

Appearance	Appearance descriptions are general information and not specific to any lot.
Physical state	Solid.
Form	Powder.
Color	White. Off-white.
Odor	Odorless.
Odor threshold	Not available.
pH	Not available.
Melting point	170 - 184 °C (338 - 363.2 °F)

Material name	Cas. N	Version #	Revision date	Issue date	Page
Clopidogrel Bisulfate	120202-66-6	2	20.05.2026	06.03.2024	4/8

Initial boiling point and boiling range	Not available.
Flash point	Not available.
Evaporation rate	Not available.
Flammability (solid, gas)	Not available.

Upper/lower flammability or explosive limits

Flammability limit - lower (%)	Not available.
Flammability limit - upper (%)	Not available.
Explosive limit - lower (%)	Not available.
Explosive limit - upper (%)	Not available.
Vapor pressure	Not available.
Vapor density	Not available.
Relative density	Not available.

Solubility(ies)

Solubility (water)	Neutral pH: Practically insoluble. Acetone: Slightly soluble.
Solubility (other)	Ethyl ether: Practically insoluble. Methanol: Freely soluble. Methylene chloride: Sparingly soluble.
Auto-ignition temperature	Not available.
Decomposition temperature	Not available.
Viscosity	Not available.

9.3 Other information

Chemical family	Thienopyridine.
Molecular formula	C ₁₆ H ₁₆ ClNO ₂ S H ₂ SO ₄
Molecular weight	419.9

10.Stability and Reactivity

Reactivity	The product is stable and non-reactive under normal conditions of use, storage and transport.
Chemical stability	Material is stable under normal conditions.
Possibility of hazardous reactions	No dangerous reaction known under conditions of normal use.
Conditions to avoid	Contact with incompatible materials.
Incompatible materials	Oxidizing agents. Acids.
Hazardous decomposition products	Irritating and/or toxic fumes or gases. Emits toxic fumes under fire conditions. NO _x . SO _x . Cl ⁻ .

11.Toxicological Information

Information on likely routes of exposure

Inhalation	Knowledge about health hazard is incomplete.
Skin contact	Causes severe skin burns.
Eye contact	Causes serious eye damage.

Material name	Cas. N	Version #	Revision date	Issue date	Page
Clopidogrel Bisulfate	120202-66-6	2	20.05.2026	06.03.2024	5/8

Ingestion	Harmful if swallowed. Based on information from therapeutic use, this material may cause: Platelet inhibition. Chest pain. Pinpoint red spots on skin. Flu-like symptoms. Cough. Shortness of breath. Swelling of feet or lower legs. Nosebleeds. Bleeding or bruising. Joint pain. Fainting. Frequent, painful, or difficult urination. Heartburn. Gastrointestinal disturbances. Back pain. Dizziness. Headache. Skin rash. Itching. Anxiety. Depression. Tiredness. Weakness. Numbness or tingling of skin. Insomnia. Leg cramps.
Symptoms related to the physical, chemical, and toxicological characteristics	

Information on toxicological effects

Acute toxicity

Product	Species	Test Results
Clopidogrel Bisulfate (CAS 120202-66-6)		
Oral LD50	Rat	1914 mg/kg
Skin corrosion/irritation	Causes severe skin burns and eye damage.	
Serious eye damage/eye irritation	Causes serious eye damage.	

Respiratory or skin sensitization

Respiratory sensitization	Knowledge about health hazard is incomplete.
---------------------------	--

Skin sensitization	Knowledge about health hazard is incomplete.
--------------------	--

Germ cell mutagenicity	Knowledge about mutagenicity is incomplete.
	Mutagenicity, In vitro Ames test Result: Negative.
	Mutagenicity, In vitro DNA repair test in rat hepatocytes Result: Negative.
	Mutagenicity, In vitro gene mutation test in Chinese hamster fibroblasts Result: Negative.
Mutagenicity	Mutagenicity, In vitro metaphase chromosome analysis of human lymphocytes Result: Negative.
	Mutagenicity, In vivo micronucleus test in mice Result: Negative.

Carcinogenicity	Based on available data, the classification criteria are not met.
-----------------	---

IARC Monographs. Overall Evaluation of Carcinogenicity	Not listed.
--	-------------

OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)	Not regulated.
--	----------------

Reproductive toxicity	Based on available data, the classification criteria are not met.
	300 mg/kg/day Reproductivity Result: Negative: No effect on fertility and no adverse effect on the fetuses. Species: Rabbit
Reproductivity	500 mg/kg/day Reproductivity Result: Negative: No effect on fertility and no adverse effect on the fetuses. Species: Rat

Specific target organ toxicity - single exposure	Causes damage to organs (Blood).
--	----------------------------------

Material name	Cas. N	Version #	Revision date	Issue date	Page
Clopidogrel Bisulfate	120202-66-6	2	20.05.2026	06.03.2024	6/8

Specific target organ toxicity - repeated exposure	Knowledge about health hazard is incomplete.
Aspiration hazard	Based on available data, the classification criteria are not met.
Further information	Pharmacologically active material. Occupational exposure may cause physiological effects.

12. Ecological information

Ecotoxicity	The product is not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment.
Persistence and degradability	No data is available on the degradability of this substance.
Bioaccumulative potential	
Octanol/water partition coefficient log K _{ow}	No data available.
Mobility in soil	No data available.
Other adverse effects	No other adverse environmental effects (e.g. ozone depletion, photochemical ozone creation potential, endocrine disruption, global warming potential) are expected from this component.

13. Disposal considerations

Disposal instructions	Dispose in accordance with all applicable regulations. Under RCRA, it is the responsibility of the user of the product to determine, at the time of disposal, whether the product meets RCRA criteria for hazardous waste.
Local disposal regulations	Dispose in accordance with all applicable regulations.
Hazardous waste code	The waste code should be assigned in discussion between the user, the producer and the waste disposal company.
Waste from residues / unused products	Empty containers or liners may retain some product residues. This material and its container must be disposed of in a safe manner.
Contaminated packaging	Since emptied containers may retain product residue, follow label warnings even after container is emptied.

14. Transport information

DOT	UN3261
IATA	UN3261
General information	It is the shipper's responsibility to determine the correct transport classification at the time of shipment.

15. Regulatory information

Safety, health and environmental regulations/legislation specific for the substance

This Safety Data Sheet has been prepared in accordance with applicable GHS principles for hazard communication.

Chemical safety assessment

A chemical safety assessment has not been carried out for this substance.

16. Other information, including date of preparation or last revision

Material name	Cas. N	Version #	Revision date	Issue date	Page
Clopidogrel Bisulfate	120202-66-6	2	20.05.2026	06.03.2024	7/8

Disclaimer

Reflabpharm CRMs are sold for chemical test and assay purposes only, and NOT for human consumption. The information contained herein applies solely to the chemical substance when used as a Reference Standard and does not necessarily relate to any other use of the substance described, (i.e. at different concentrations, in drug dosage forms, or in bulk quantities). CRMs are intended for use by persons having technical skills and at their discretion and risk. This information has been developed by Reflabpharm staff from sources considered reliable but has not been independently verified by Reflabpharm. Therefore, the Reflabpharm cannot guarantee the accuracy of the information in these sources nor should the statements contained herein be considered an official expression. NO REPRESENTATION OR WARRANTY, EXPRESS OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE is made concerning the information contained herein.

Material name	Cas. N	Version #	Revision date	Issue date	Page
Clopidogrel Bisulfate	120202-66-6	2	20.05.2026	06.03.2024	8/8